





# Tende VAV-IW,

Infant Warmer is a comprehensive infant care system that also supports the clinicians with its advanced algorithm and utmost automation feature to begin to act instantly when necessary.

## Color TFT Touch Screen Display

Tende VAV-IW, Infant Radiant Warmer, provides users with accurate, effective and fast operation and easy monitoring of patient data with 12.1 inch color TFT touch screen with graphics, animations and detailed explanations.

With the 12.1 inch screen, users can easily observe intense data such as the Masimo Rainbow SET pulse oximeter (SpO2) and non-invasive blood pressure measurement (NIBP). With the Easy Mute feature, the clinician's alarms that occur during the procedure are silenced with just a hand gesture without touch.





### Lighting

Clinicians have the opportunity to observe the baby by using examination lighting and/or spot lighting according to their needs.

## Electronic Mattress Tilt (Trendelenburg)

The user can easily adjust the inclination of the bed electronically without manpower on the color touch screen. This process, which is done silently and without vibration, does not disturb the baby. The fact that the bed can be brought to the horizontal position with a single button enables the clinician to reach the baby without wasting time.



## **Membrane Control Panel**

The clinician can easily give commands and monitor data with membrane keys on the control panel suitable for usage habits.



## $360^{\circ}$ Rotatable Bed

The 360° rotatable mattress and the heater unit that can be rotated in both directions allow surgical operations and x-rays imaging to be performed easily.

#### Large Drawers

Tende VAV-IW has a large drawer with an organizer that can be accessed from the right and left sides of the infant warmer.



## CONFIGURABLE STRUCTURE

Tende VAV-IW, Infant Warmer can be configurable from the basic model to advanced features like resuscitator, air/oxygen blender, suction module, height adjustment, phototherapy, pulse oximeter, NIBP, digital scale.

## **Resuscitator (Optional)**

The newborn can have consistent PIP and PEEP values for efficient oxygenation. There are two more option Air/Oxygen Blender and Suction Module

## Integrated Wireless Digital Scale (Optional)

Tende VAV-IW high sensitive weighing scale was designed considering the safety and comfort of the infant. Just lifting the newborn for a few seconds is enough to tare, then the system performs weighing operations when the infant lowered. Operation stages can be guided by animation from the screen with sound.

## NIBP (Optional)

NIBP (Non-Invasive Blood Pressure Measurement); Blood pressure and pulse values for a infant can be measured with determined time intervals. Clinicians can determine set values and activate audible and visual alarms.

## Pulse-Oximeter (Optional)

Tende VAV-IW uses the most advanced Masimo Rainbow SET SpO2 technology. The advantage of this set is that it is also possible to measure perfusion index (PI) and pleth variability index (PVI), total hemoglobin (SpHb) methemoglobin (SpMet) carboxyhemoglobin (SpCO) values without taking blood samples from the newborn (Non-Invasive).

## Electronic Height Adjustment (Optional)

Clinicians can adjust the height of the infant warmer for a comfortable working height. Lifting and lowering operations can be made silently without vibration and without disturbing the newborn by means of the motor which is specially produced for Tende VAV-IW. This operation can be implemented by foot pedals from both sides of the infant warmer.



#### **TECHNICAL SPECIFICATIONS**



ENERAL SPECIFICATIONS	
isplay	12.1 inch Color TFT Touch Screen Display
izes	Height (Standard): 182 cm, Width: 63 cm, Depth: 115 cm
lectronic Height Adjustment	Mattress Height; Min. 90 cm (± 1cm) Max. 110 cm (± 1cm)
Veight Vheels	Standard;120 Kg, (Without options) 4 Dual Wheeled Casters with All Friction Brake (Diameter: 125 mm)
Aonitorization Backup	Monitorization After Power Failure (Min. 20 Minutes)
larm Sound Level Adjustment	Yes
ower Requirement	220-240V ± 10%, 50/60 Hz
PERATING CONDITIONS	
emperature	15°C to 40°C
lumidity	0% to 90% (Noncondensing Relative Humidity)
NODES	Add stable O Ma Alexan
rewarm Mode	Adjustable & No Alarms Automated
aby Mode Ianuel Mode	User Controlled
ASSINET SPECIFICATIONS	User controlled
fattress	Biocompatible, Cytotoxicity and Skin Irritation Free, Fire Retarded,
	Water Proof & Washable (ISO 10993-10 Sensitivity Reported, ISO 10993-5)
Nattress Size	Length; 72 cm x Width; 55 cm x Height; 3 cm
1attress Tilt (Electronic Control)	+ 15° to -12°, Continuously Variable
otatable Bed	360°
rommets	4 Pcs
rawer	Large drawer with organizer and file holder (Two way opening)
ube Holders ( Polo & Monitor Table	Yes (2pcs)(Optional)
/ Pole & Monitor Table EATER & TEMPERATURE	Yes
uto Heater	Yes
emperature Gauge	15°C - 42°C
kin Temperature Control	34°C to 38°C with 0.1°C Resolution )(>37°C with confirm)
eramic Heater Life Time	> 20,000 Hours
THER FEATURES	·
mer	Patented; Timer & Countdown (Timer also APGAR, CPR)
xamination Lamp (LED)	Yes
pot Examination Lamp (LED)	Yes
PGAR Timer	Yes
eater Head Rotation	Yes 180°
ED PHOTOTHERAPY (OPTIONAL)	
ght Source	Blue LED (20 pcs, Super Bright LED)
/avelength ED Life Time	460nm (± 2%) > 70,000 Hours (Check the user manual for technical details)
ight Intensity (Standard)	Adjustable
ESUSCITATION MODULE (OPTIONAL)	Aujustable
eak Inspiration Pressure (PIP)	Yes
ositive-end Expiratory Pressure(PEEP)	Yes
low Capacity	Yes
irway Pressure Manometer	Yes
IR/O2 BLENDER	Yes, Optional
ange	21-100 %O2
ccuracy	± 3% full scale
lowmeter Range	0-15 lpm
UCTION	Yes, Optional
/aste Jar Capacity	500 ml
ITEGRATED DIGITAL SCALE (OPTIONAL)	100  ar to  10  kg(10,000  ar)
/eight Range	100 gr to 10 kg (10,000 gr)
/eight Accuracy /eight Display Resolution	± 5gr 1 gr
DO2 & PULSE OPTION (Masimo Rainbow SET)	± δ'
002 Display Range	0-100% SpO2 with 1% Resolution
ulse Rate Display Range	0 to 240 BPM with 1 BPM Resolution
, Arterial Pulse Signal Strength	
	Perfusion Index 0 to 20 Values
	0 to 100 Values
VI, Pleth Variability Index oMet, Methemoglobin Saturation	0 to 100 Values 0% to 100%
/I, Pleth Variability Index DMet, Methemoglobin Saturation DCO, Carboxyhemoglobin Saturation	0 to 100 Values 0% to 100% 0% to 100%
VI, Pleth Variability Index DMet, Methemoglobin Saturation DCO, Carboxyhemoglobin Saturation DHb, Total Hemoglobin Concentration in Arterial Blood	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L
VI, Pleth Variability Index DMet, Methemoglobin Saturation DCO, Carboxyhemoglobin Saturation DHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index	0 to 100 Values 0% to 100% 0% to 100%
VI, Pleth Variability Index oMet, Methemoglobin Saturation oCO, Carboxyhemoglobin Saturation oHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric)	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values
/I, Pleth Variability Index oMet, Methemoglobin Saturation oCO, Carboxyhemoglobin Saturation oHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric) ood Pressure Range	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg
/I, Pleth Variability Index DMet, Methemoglobin Saturation DCO, Carboxyhemoglobin Saturation DHD, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index <b>IBP OPTION (Method of Measurement; Oscillometric)</b> ood Pressure Range uff Pressure	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C)
VI, Pleth Variability Index DMet, Methemoglobin Saturation DCO, Carboxyhemoglobin Saturation DHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric) ood Pressure Range uff Pressure uff Deflate Rate	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume
/I, Pleth Variability Index DMet, Methemoglobin Saturation DCO, Carboxyhemoglobin Saturation DHD, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric) ood Pressure Range Uff Pressure uff Deflate Rate itial Inflation Pressure	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg
/I, Pleth Variability Index Met, Methemoglobin Saturation OCO, Carboxyhemoglobin Saturation Hb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric) ood Pressure Range Iff Pressure Iff Deflate Rate itial Inflation Pressure Iff-Adjusting Inflation Pressure	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume
VI, Pleth Variability Index bMet, Methemoglobin Saturation bCO, Carboxyhemoglobin Saturation bHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index <b>IBP OPTION (Method of Measurement; Oscillometric)</b> lood Pressure Range uff Deflate Range uff Deflate Rate itial Inflation Pressure elf-Adjusting Inflation Pressure inical Accuracy (Meets accuracy requirements of)	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg
VI, Pleth Variability Index bMet, Methemoglobin Saturation bCO, Carboxyhemoglobin Saturation bHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric) lood Pressure Range uff Deflate Rate itial Inflation Pressure elf-Adjusting Inflation Pressure inical Accuracy (Meets accuracy requirements of) cartup Initialization Period	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013
VI, Pleth Variability Index bMet, Methemoglobin Saturation bCO, Carboxyhemoglobin Saturation bHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric) ood Pressure Range uff Pressure uff Deflate Rate itial Inflation Pressure elf-Adjusting Inflation Pressure inical Accuracy (Meets accuracy requirements of) cartup Initialization Period ong Term Interval Periods REND PARAMETERS	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes
VI, Pleth Variability Index bMet, Methemoglobin Saturation bCO, Carboxyhemoglobin Saturation bCD, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IBP OPTION (Method of Measurement; Oscillometric) ood Pressure Range uff Pressure uff Deflate Rate itial Inflation Pressure lef-Adjusting Inflation Pressure lef-Adjusting Inflation Pressure inical Accuracy (Meets accuracy requirements of) cartup Initialization Periods REND PARAMETERS 8, 24, 168 hours (7 Days) Trend Parameters for	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds
VI, Pleth Variability Index bMet, Methemoglobin Saturation bCO, Carboxyhemoglobin Saturation bHD, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index <b>IBP OPTION (Method of Measurement; Oscillometric)</b> lood Pressure Range uff Pressure uff Deflate Rate itial Inflation Pressure elf-Adjusting Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) cartup Initialization Periods <b>REND PARAMETERS</b> . 8, 24, 168 hours (7 Days) Trend Parameters for <b>LARMS (Audible &amp; Visual)</b>	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes Skin Temperature , Weight, Heater Power, NIBP
VI, Pleth Variability Index oMet, Methemoglobin Saturation oCO, Carboxyhemoglobin Saturation oHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index <b>IBP OPTION (Method of Measurement; Oscillometric)</b> lood Pressure Range uff Deflate Rate uff Deflate Rate itial Inflation Pressure elf-Adjusting Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) cartup Initialization Period ong Term Interval Periods <b>REND PARAMETERS</b> . 8, 24, 168 hours (7 Days) Trend Parameters for <b>LARMS (Audible &amp; Visual)</b> aby Check Alarm	0 to 100 Values         0% to 100%         0% to 100%         0% to 25% g/dL, 0% to 16% mmol/L         0.00 to 1.00 Values         Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg         Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C)         Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume         90 mmHg (Default), Variable from 60 to 140 mmHg         Previous Systolic + 20 mmHg         ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013         6 Seconds         1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes         Skin Temperature , Weight, Heater Power, NIBP         Yes
VI, Pleth Variability Index pMet, Methemoglobin Saturation pCO, Carboxyhemoglobin Saturation pHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IIBP OPTION (Method of Measurement; Oscillometric) lood Pressure Range uff Pressure uff Deflate Rate nitial Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) tartup Initialization Period ong Term Interval Periods REND PARAMETERS , 8, 24, 168 hours (7 Days) Trend Parameters for LARMS (Audible & Visual) aby Check Alarm kin (In Baby Mode)	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes Skin Temperature , Weight, Heater Power, NIBP Yes High/Low Temp., Over Temp., Skin Probe Read Fault, Skin Probe Disconnected
VI, Pleth Variability Index pMet, Methemoglobin Saturation pCO, Carboxyhemoglobin Saturation pHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IIBP OPTION (Method of Measurement; Oscillometric) lood Pressure Range uff Deflate Rate nitial Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) tartup Initialization Period ong Term Interval Periods REND PARAMETERS , 8, 24, 168 hours (7 Days) Trend Parameters for LARMS (Audible & Visual) aby Check Alarm kin (In Baby Mode) pO2, Pulse Rate	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes Skin Temperature , Weight, Heater Power, NIBP Yes High/Low Temp., Over Temp., Skin Probe Read Fault, Skin Probe Disconnected High/Low Values, Related Security Alarms
VI, Pleth Variability Index pMet, Methemoglobin Saturation pCO, Carboxyhemoglobin Saturation pHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index <b>IBP OPTION (Method of Measurement; Oscillometric)</b> lood Pressure Range uff Deflate Rate nitial Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) tartup Initialization Period ong Term Interval Periods <b>REND PARAMETERS</b> , 8, 24, 168 hours (7 Days) Trend Parameters for <b>LARMS (Audible &amp; Visual)</b> aby Check Alarm kin (In Baby Mode) pO2, Pulse Rate	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes Skin Temperature , Weight, Heater Power, NIBP Yes High/Low Temp., Over Temp., Skin Probe Read Fault, Skin Probe Disconnected
VI, Pleth Variability Index pMet, Methemoglobin Saturation pCO, Carboxyhemoglobin Saturation pHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IIBP OPTION (Method of Measurement; Oscillometric) lood Pressure Range uff Deflate Rate nitial Inflation Pressure elf-Adjusting Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) tartup Initialization Period ong Term Interval Periods REND PARAMETERS , 8, 24, 168 hours (7 Days) Trend Parameters for LARMS (Audible & Visual) aby Check Alarm kin (In Baby Mode) pO2, Pulse Rate thers	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes Skin Temperature , Weight, Heater Power, NIBP Yes High/Low Temp., Over Temp., Skin Probe Read Fault, Skin Probe Disconnected High/Low Values, Related Security Alarms Heater Temp. Fault, Power, Low Battery, Excess Weight, System Failure etc.
VI, Pleth Variability Index pMet, Methemoglobin Saturation pCO, Carboxyhemoglobin Saturation pHb, Total Hemoglobin Concentration in Arterial Blood PRI, Oxygen Reserve Index IIBP OPTION (Method of Measurement; Oscillometric) lood Pressure Range uff Deflate Rate nitial Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) tartup Initialization Period ong Term Interval Periods REND PARAMETERS , 8, 24, 168 hours (7 Days) Trend Parameters for LARMS (Audible & Visual) aby Check Alarm kin (In Baby Mode) pO2, Pulse Rate thers 1984 MDD 93/42/EEC ANNEX 2 / SEC 3	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes Skin Temperature , Weight, Heater Power, NIBP Yes High/Low Temp., Over Temp., Skin Probe Read Fault, Skin Probe Disconnected High/Low Values, Related Security Alarms
VI, Pleth Variability Index pMet, Methemoglobin Saturation pCO, Carboxyhemoglobin Saturation pHb, Total Hemoglobin Concentration in Arterial Blood Ri, Oxygen Reserve Index IIBP OPTION (Method of Measurement; Oscillometric) lood Pressure Range uff Deflate Rate nitial Inflation Pressure elf-Adjusting Inflation Pressure elf-Adjusting Inflation Pressure linical Accuracy (Meets accuracy requirements of) tartup Initialization Period ong Term Interval Periods REND PARAMETERS , 8, 24, 168 hours (7 Days) Trend Parameters for LARMS (Audible & Visual) aby Check Alarm kin (In Baby Mode) pO2, Pulse Rate thers	0 to 100 Values 0% to 100% 0% to 100% 0% to 25% g/dL, 0% to 16% mmol/L 0.00 to 1.00 Values Systolic 40-130 mmHg, MAP 26-110 mmHg, Diastolic 20-100mmHg Measurement Range; 0-300 mmHg, Accuracy; ± 3 mmHg (0-50°C) Deflation Step Size Varies with Heart Rate, Cuff Pressure and Cuff Volume 90 mmHg (Default), Variable from 60 to 140 mmHg Previous Systolic + 20 mmHg ANSI/AAMI SP10:2002(R)2008, EN1060-4:2004 and ISO 81060-2:2009/2013 6 Seconds 1, 2, 3, 4, 5, 10, 15, 30, 60, 90, 120, 240 Minutes Skin Temperature , Weight, Heater Power, NIBP Yes High/Low Temp., Over Temp., Skin Probe Read Fault, Skin Probe Disconnected High/Low Values, Related Security Alarms Heater Temp. Fault, Power, Low Battery, Excess Weight, System Failure etc.